

Medi One

Distributeur exclusif France



Single Use Flexible Video Uretero-
Renscope

Specifications technique

Optical System	Pixels: CMOS 400x400 soit 160K
	Champ de vision : 120°
	Plage visible: 3-50mm
	Eclairage: Fibre optique LED
Insertion Tube	Diametre externe max : 9Fr
	longueur du cordon : 630mm
Instrument Channel	Diametre du canal de travail :3.6Fr
Bending Section	Flexion : 280° dans les deux sens
Function	Photo
	Video
	Balance des blancs
	Enregistrement : 16Go disponible
Connections	AV
	HDMI
	VGA
	DVI
	SDI

AN ENDOSCOPE SOLUTION COMPANY

No Compromise
No Repairs
No Sterilization
Always Available
Cost Effective



Shenzhen HugeMed Medical Technical Development Co.,Ltd

Add: 8F,Block A, Tempus Building,1st Qingshuihe Road,
Luohu District Shenzhen City, Guangdong, China
Tel: +86-755-2227 5866
+86-755-8226 7820
Fax: +86-755-2227 5833
Email: sale@hugemed.net


www.hugemed.net



EndoView HU30

Single-Use Uretero-Renoscope

Hugemed all in one system, always be ready for portable image processor



285°
Excellent Deflection

For easier navigation of lower pole

- Challenging Anatomy
- Passive Bending



160K High Resolution



Dual Buttons
For snap / Video capture



Self Locker
For reducing surgeon fatigue



True 1: 1
Torque ratio

Flexible Sheath
3.6Fr working channel