

Medi One

Distributeur exclusif France



HU30S Single use Uretero-Renoscope

Technical Specifications

| | |
|--------------------|-------------------------------------|
| Optical System | Pixels: CMOS 400x400 soit 160K |
| | Champ de vision : 120° |
| | Plage visible: 3-50mm |
| | Eclairage: Fibre optique |
| Insertion Tube | Diametre externe max : 9Fr |
| | longueur du cordon : 630mm |
| Instrument Channel | Diametre du canal de travail :3.6Fr |
| Bending Section | Flexion : 280° dans les deux sens |

Smaller and Safer

Shenzhen HugeMed
Medical Technical Development Co.,Ltd

Add: 8F,Block A, Tempus Building,1st Qingshuihe Road,
Luohu District Shenzhen City, Guangdong, China
Tel: +86-755-2227 5866
+86-755-8226 7820
Fax: +86-755-2227 5833
Email: sale@hugemed.net

www.hugemed.net

AN ENDOSCOPE SOLUTION COMPANY

No Compromise
No Repairs
No Sterilization
Always Available
Cost Effective




7.5 Fr


HU30S SINGLE-USE URETERO- RENSCOPE

RIRS


- Reliable irrigation flow
- Increased stone free rate
- Reduced intrarenal pressure
- Safer ureter access



- Same excellent visualization
- Same excellent deflection
- Same working channel dimension



285°
Excellent Deflection
For easier navigation
of lower pole



The HU30S offers excellent visualization even in the most confined spaces.

